

PachPen[®]

Handheld Pachymeter

Item #: 24-5100

Features

- Accurate – 65 MHz sampling probe ensures precise measurements
- Gentle Touch – 2.5mm diameter probe helps to ensure patient comfort
- Long Lasting Battery – Long lasting lithium battery lasts for approximately 15,000 readings
- IOP Correction Calculator – Converts IOP measurements in seconds
- 9 Readings – With a calculated average

Specifications

ACCURACY	+/- 5 Microns
WEIGHT	3.0 oz.
DIMENSIONS	7 ¼" x 1 ¼"
PROBE DIMENSION	2.5 mm Diameter
MEASUREMENT RANGE	300–999 Microns
POWER SOURCE	3.6 V Lithium Battery



Accessories Included

- (1) Hard shell carrying case
- (2) XENO 3.6V Lithium Battery
- (1) User's Guide
- (1) Screwdriver
- (1) Lanyard
- (24) Alcohol prep pads

Please contact your Keeler sales representative for additional information.